

Ultra-sensitive detection and quantification of HPV DNA in the plasma of patients with oropharyngeal squamous cell carcinoma (OPSCC) enrolled in the OPTIMA 2 treatment de-escalation trial

Hillary Sloane¹, Evgeny Izumchenko², Austin Mattox³, Rifat Hasina², Aashay Patel⁴, Frederick S. Jones¹, Hannah Quinn¹, Kimberly Schmitt¹, Jennifer Preston¹, Kirsten Keyser⁵, Johannes Fredebohm⁵, Frank Holtrup⁵, Alexander T. Pearson², Everett E. Vokes², Tanguy Y. Seiwert³, Daniel L. Edelstein¹, Ari Rosenberg², Nishant Agrawal²

¹Sysmex Inostics, Inc., Baltimore, MD, USA; ²University of Chicago, Chicago, IL, USA; ³Johns Hopkins University, Baltimore, MD USA; ⁴Icahn School of Medicine at Mount Sinai, New York, NY, USA; ⁵Sysmex Inostics GmbH, Hamburg, Germany



Introduction

- Human papillomavirus (HPV) infection is a primary factor driving the increasing incidence of OPSCC.
- As patients with HPV+ OPSCC show significantly improved treatment response and prognosis, there is an urgent need to de-escalate treatment of HPV+ OPSCC that optimizes oncologic control while minimizing treatment-related toxicity.
- Cell-free HPV DNA (cfHPV-DNA)** from plasma represents a promising noninvasive surrogate of disease burden in these patients.
- To enable cfHPV-DNA analysis as a strategy to monitor therapy response and guide treatment de-escalation, we developed a highly sensitive assay for HPV16/18 detection and quantification in plasma, based on SafeSEQ next-generation sequencing (NGS) technology.

Methods

- Longitudinal plasma samples were collected from patients with locoregionally HPV+ OPSCC treated on OPTIMA 2 (NCT03107182), a de-escalation protocol of induction chemioimmunotherapy followed by risk/response stratified de-escalated locoregional therapy.
- Neck CT or MRI was obtained for all patients at baseline and following induction chemioimmunotherapy; radiographic response to induction therapy was assessed per RECIST 1.1 criteria.
- cfHPV-DNA was quantified in plasma samples (2mL plasma per time-point) collected at baseline and at 6-9 weeks after beginning induction therapy using the SafeSEQ cfHPV-DNA Test (Sysmex Inostics).





Figure 1. OPTIMA 2 study design. TORS: Transoral Robotic Surgery; TFHX: paclitaxel, 5-fluorouracil, and hydroxyurea, with twice daily radiation in week-on week-off CRT platform

Conclusions

- The SafeSEQ cfHPV-DNA Test exhibits robust quantitative detection of HPV across a broad range, enabling high-resolution molecular monitoring for HPV+ OPSCC patients.
- Prospective studies are underway to further evaluate the kinetics of cfHPV-DNA as a predictor of response to therapy in order to more precisely guide the management of patients with HPV+ OPSCC.

Results


SafeSEQ cfHPV-DNA Test: Analytical performance



- Accurate quantification of plasma HPV 16 and 18 levels over 5 orders of magnitude.
- Low quantitative variability (<30% CV observed above 5 copies).
- Low level of background signal (<0.04 copies per sample across 20 healthy donor samples), indicating high specificity.

Figure 2. Quantification of HPV 16 and 18 in contrived samples. Dilution series of contrived samples tested in replicate at 12 tiers ranging from 0.3 to 50,000 copies

cfHPV-DNA detected in patient plasma samples collected at baseline



- In 32 (82%) of 39 patient plasma samples collected at baseline, cfHPV-DNA was detected at levels ranging from 2.5 to >59,000 copies.
- cfHPV-DNA was not detected above the assay cut-off (1.5 copies) in 7 of the 39 baseline samples tested.

Figure 3. cfHPV-DNA detected in 32 patient plasma samples collected at baseline.

Longitudinal cfHPV-DNA analysis

- Plasma samples collected at 2 serial time-points (baseline and 6-9 weeks after beginning induction therapy) were available for analysis in 31 patients. cfHPV-DNA was detected at baseline in 25 patients.
- All 25 patients with baseline cfHPV-DNA showed a decrease in cfHPV-DNA level at follow-up, with complete clearance observed in 20/25 (80%) of patients; this is consistent with tumor response (shrinkage) to induction therapy observed in 24/25 (96%) patients.
- Of 5 patients with persistent cfHPV-DNA detected post-induction therapy (indicated by yellow highlighted rows in Table 1): 1 patient progressed on induction therapy, 1 patient demonstrated subsequent recurrence and death, and 1 patient demonstrated concern for distant metastasis followed by death.




Figure 4. Longitudinal cfHPV-DNA analysis. Trend lines are colored according to radiographic response to induction therapy. Dashed lines indicate patients with persistent cfHPV-DNA at follow-up; residual cfHPV-DNA data points are indicated with a ♦

Patient ID	Clinical stage	Clinical risk	Baseline cHPV16	Follow-up sample		Post-induction response (% tumor shrinkage)	Clinical follow-up / outcome
				cHPV16	Weeks from baseline		
UC005	T3N0	LR	4,622.75	0.00	9	100	NER
UC017	T2N2b	LR	302.63	0.42	8	100	NER
UC023	T4N2c	HR	180.45	0.00	9	100	NER
UC064	T2N2b	HR	71.73	0.00	8	100	NER
UC030	T2N2b	LR	30.09	0.00	8	100	NER
UC019	T3-4N2b	HR	0.05	0.00	8	100	NER
UC011	T3N1	HR	76.57	1.83	9	87.7	NER
UC014	T4N2c	HR	6,173.74	1.97	9	81.6	Concern for lung met (deceased)
UC016	T4N3	HR	9,740.48	0.00	9	74.6	NER
UC027	T2N2c	HR	401.31	0.68	8	73.3	NER
UC008	T1N2	HR	694.18	0.00	7	68.8	NER
UC024	T1N2a	LR	51.71	0.33	8	65.7	NER
UC013	T2N2b	LR	5,705.97	0.00	8	64.5	NER
UC067	T1N2b	LR	365.26	0.00	8	63.9	NER
UC068	T2N1	HR	12.88	0.00	9	63.0	NER
UC028	T1N2b	HR	12.86	0.00	9	60.7	NER
UC062	T2N2b	LR	16.72	0.00	8	55.3	NER
UC018	T2N1	LR	38.44	0.00	9	55.1	NER
UC066	T1N2b	LR	0.38	0.00	8	54.6	NER
UC015	T2N2b	LR	1.09	0.00	8	53.3	NER
UC009	T1N2a	LR	985.17	0.00	8	51.9	NER
UC022	T2N2b	HR	0.10	0.00	8	50.7	NER
UC070	T1N2b	HR	946.06	0.00	8	48.0	NER
UC020	T2N2b	LR	54.77	6.50	8	46.9	Recurrence (deceased)
UC003	T1N2	LR	37,855.87	0.00	8	45.5	NER
UC063	T1N2b	LR	215.82	0.00	8	45.5	NER
UC001	T2N2a	LR	0.16	0.00	8	38.5	NER (deceased)
UC039	T2N2b	HR	59,050.56	0.75	7	37.5	NER
UC040	T2N2	LR	184.92	44.03	5	30.0	NER
UC007	T2N2b	HR	5.32	0.00	8	28.6	NER
UC073	T4N2b	HR	1,084.91	19.94	6	2% GROWTH	NER

Table 1. Summary of patients included in longitudinal cfHPV-DNA analysis. HPV values in grey italics are below assay cut-off; LR: Low risk; HR: High risk; NER: No evidence of recurrence.